

TACTICAL ASSAULT BENCH FOR THE BELL 407 *SERIES II* FOLDING ALWiP



Product Information – *TACTICAL ASSAULT BENCH ATTACHMENT*



Tactical Assault Bench (TAB) on the BELL 407 *Series II* COMPOSITE BASED Folding platform in aircraft shell

Cantine Armament Inc, proprietary rights are included in the information disclosed herein. Recipient by accepting this document agrees that neither this document nor the information disclosed herein nor any part thereof shall be reproduced or transferred to other documents or used or disclosed to others for manufacturing or any other purpose except as specifically authorized in writing by Cantine Armament Inc.



TACTICAL ASSAULT BENCH FOR THE BELL 407 *SERIES II* FOLDING ALWiP



Product Information – *TACTICAL ASSAULT BENCH ATTACHMENT*


The new Tactical Assault Bench (TAB) for the *Series II* Folding Advanced Lightweight Weapons Platform (ALWiP) is an elegant and dependable solution for rapid operator insertion/extraction. The TAB is constructed from the same composite materials as the *Series II* ALWiP to maintain strength and unity within the system. The lightweight modular system is designed to be quickly installed, removed, or reconfigured to support rapid deployments and flexible mission profiles.

The Tactical Assault Bench has the capability to transport up to three fully-equipped operators per bench and can be configured with an optional pintle mount for adding various weaponry to include the M-240, M-249, Mk 48, and M-134 mini gun. Our standard cargo rail system is incorporated for the attachment of operator safety lines, casualty litters, etc. The TAB is designed with a tactile lip and scalloped profile to provide intuitive positioning feedback to the operator without losing focus on the mission. The system features an optional removable footrest and a forward wind deflector to improve operator comfort at high speeds.

The modular design allows for the addition of each bench to be attached to the center portion of the platform using the existing hardware in less than two minutes with no additional tools. Alternately, the TAB may be quickly removed and the aircraft doors secured with no additional reconfiguring required. These capabilities allow for rapid and flexible deployment of the helicopter as mission requirements change.

Features:

- Modular lightweight design allows for rapid installation
- Standard cargo rail system for additional mission configurations
- Three fully-equipped operator capacity per bench
- Optional removable footrest and forward wind deflectors
- One litter patient per bench



Tactical Assault Bench (TAB) on the BELL 407 *Series II* ALWiP

Cantine Armament Inc, proprietary rights are included in the information disclosed herein. Recipient by accepting this document agrees that neither this document nor the information disclosed herein nor any part thereof shall be reproduced or transferred to other documents or used or disclosed to others for manufacturing or any other purpose except as specifically authorized in writing by Cantine Armament Inc.